



# Recombinant Arabidopsis thaliana Alpha-ketoglutarate-dependent dioxygenase alkB (At1g11780)

<b>Product Code</b>	CSB-EP874366DOA
<b>Abbreviation</b>	At1g11780
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9SA98
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MYESANVSDD ADRTAFRRAE KKYKLYYEQD SKFSRKKKLP KPIDLSELLD FNLISQNFNN DGVLPDGIRV SKVDSSPVFC IDNRPGFYFI PDALSLKEQC KWIKESLTSF PQPPNRTNHN AIYGPIDDLF DSAKENKVLV QDDLTTNNKWK FYEEVDIEKA TRSSCKSVSA SVLLRKLKRW TLGLQFDWSK RNYDVSLPHN NIPDALCQLA KTHAAIAMPD GEEFRPEGAI VNYFGIGDTL GGHLDDMEAD WSKPIVSMSL GCKAIFLLGG KSKDDPPHAM YLRSGDVVLM AGEARECFHG IPRIFTGEEN ADIGALESEL SHESGHFFAE YIKTSRININ IRQVF
<b>Source</b>	E.coli
<b>Target Names</b>	At1g11780
<b>Protein Names</b>	Recommended name: Alpha-ketoglutarate-dependent dioxygenase alkB EC= 1.14.11.- Alternative name(s): Alkylated DNA repair protein alkB homolog
<b>Expression Region</b>	1-345
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.