



Recombinant Arabidopsis thaliana Alpha-ketoglutarate-dependent dioxygenase alkB (At1g11780)

Product Code	CSB-YP874366DOA
Abbreviation	At1g11780
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SA98
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MYESANVSDD ADRTAFRRAE KKYKLYYEQD SKFSRKKKLP KPIDLSELLD FNLISQNFNN DGVLPDGIRV SKVDSSPVFC IDNRPGFYFI PDALSLKEQC KWIKESLTSF PQPPNRTNHN AIYGPIDDLF DSAKENKVLV QDDLTTNNKWK FYEVDIEKA TRSSCKSVSA SVLLRKLKRW TLGLQFDWSK RNYDVSLPHN NIPDALCQLA KTHAAIAMPD GEEFRPEGAI VNYFGIGDTL GGHLDDMEAD WSKPIVSMSL GCKAIFLLGG KSKDDPPHAM YLRSGDVVLM AGEARECFHG IPRIFTGEEN ADIGALESEL SHESGHFFAE YIKTSRININ IRQVF
Source	Yeast
Target Names	At1g11780
Protein Names	Recommended name: Alpha-ketoglutarate-dependent dioxygenase alkB EC= 1.14.11.- Alternative name(s): Alkylated DNA repair protein alkB homolog
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.