



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At1g13200 (At1g13200)

Product Code	CSB-EP874369DOA
Abbreviation	At1g13200
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SAF4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MKDAEKREVI ASSSLQRKRN RGRRLRKRRR RNEKRVLMVP SSLPNDVLEE IFLRFVVKAL IRLKSLSKQW RSTIESRSFE ERHLTIAKKA FVDHPKVMLV GEEDPIRG TG IRPDTDIGFR LFCLESASLL SFTRLNFPQG FFWIYISES CDGLFCIHSP KSHSVYVNP ATRWLRLPP AGFQILHKF NPTEREWNVV MKSIFHLAFV KATDYKLVWL YNCDKYIVDA SSPNVGVTKC EIFDFRKNW RYLACTPSHQ IFYYQKPASA NGSVYWFTEP YNERIEVAVF DIQTETFRLL PKINPAIAGS DPHHIDMCTL DNSLCMSKRE KDTMIQDIWR LKPSEDTWEK IFSIDLVSCP SSRTEKRDQF DWSKKDRVEP ATPVAVCKNK KILLSHRYSR GLVKYDPLTK SIDFFSGHPT AYRKVIYFQS LISHL
Source	E.coli
Target Names	At1g13200
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At1g13200
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.