



Recombinant Arabidopsis thaliana F-box/kelch-repeat protein SKIP6 (SKIP6)

Product Code	CSB-MP874434DOA
Abbreviation	SKIP6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SJ04
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAATTSSGDE PPETKSPAQL IPLLSEVAL SCLARVPRCH YPILSLVSKT FRSLPTSPLL YATRALVGAT ENILYVAIRI PPESGACWFT LLHRTLNSST NSKMLVPIPS CPSPSLVGSA YVVVDSEIYV IGGIRDVPS SSVWVLCRF HTWRRVSNMR VGREFAAAGV IDGKIYVIGG CVVDNWARSI NWAEMFDIKT QTWEPVASPG MEVREKWMHA SAVMEGKVYA MADRNGVVYE PKEKKWEMPE KRLDLGWRGR ACVIENILYC YDYLKIRGY DPKERIWREL KGVESLPKFL CGATMANRGG KLTVLWEGKA GSGGSRMEI WCAEIDVGRR GEREIWGKIE WSGTVLTPN ESAIVNCLAA TV
Source	Mammalian cell
Target Names	SKIP6
Protein Names	Recommended name: F-box/kelch-repeat protein SKIP6 Alternative name(s): SKP1-interacting partner 6
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.