



Recombinant Arabidopsis thaliana Ras-related protein RABG2 (RABG2)

Product Code	CSB-MP874435DOA
Abbreviation	RABG2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SJ11
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDSLKNRTLL KVIVLGDSGV GKTSLMNQYV YKKFNKQYKA TIGADFVTKELHIDEKSVTL QIWDTAGQER FQSLGAAFYR GADCCVLVYD VNNLKSFETLNNWHITEFLKQ ANPMEPETFP FVLIGNKTDV DGGNSRVVSN KRAIEWCGSKGNIPYHETSA KEDTNIDEAF LSAHIALSN ERKQSNDIYP RGQYHDSVTD IIDPDQSRGC AC
Source	Mammalian cell
Target Names	RABG2
Protein Names	Recommended name: Ras-related protein RABG2 Short name= AtRABG2 Alternative name(s): Ras-related protein Rab77 Short name= AtRab77 Ras-related protein Rab7a Short name= AtRab7a
Expression Region	1-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.