



Recombinant Arabidopsis thaliana Ethylene-responsive transcription factor RAP2-7 (RAP2-7)

Product Code	CSB-MP874447DOA
Abbreviation	RAP2-7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SK03
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MLDLNLNADS PESTQYGGDS YLDRQTS DNS AGNRVEESGT STSSVINADG DEDSCSTRAF TLSFDILKVG SSSGGDESPA ASASVTKEFF PVSGDCGHLR DVEGSSSSRN WIDLSFDRIG DGETKLVTPV PTPAPVPAQV KKSRRGPRSR SSQYRGVTFY RRTGRWESHI WDCGKQVYLG GFDTA HAAAR AYDRAAIKFR GVDADINFTL GDYEEDMKQV QNLSKEEFVH ILRRQSTGFS RGSSKYRGVT LHKCGRWEAR MGQFLGKKAY DKAINTNGR EAVTNFEMSS YQNEINSESN NSEIDLNLGI SLSTGNAPKQ NGRLFHFP SN TYETQRGVSL RIDNEYMGKP VNTPLPYGSS DHRLYWNGAC PSYNNPAEGR ATEKRSEAE G MMSNWGWQRP GQTS AVR PQP PGPQPPLFS VAAASSGFSH FRPQPPNDNA TRGYFY PHP
Source	Mammalian cell
Target Names	RAP2-7
Protein Names	Recommended name: Ethylene-responsive transcription factor RAP2-7 Alternative name(s): Protein RELATED TO APETALA2 7 Protein TARGET OF EAT 1
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.