



Recombinant Arabidopsis thaliana Ras-related protein RABA1d (RABA1D)

Product Code	CSB-BP874475DOA
Abbreviation	RABA1D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SN35
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MAGYRADDDY DYLFKVVLLIG DSGVGKSNLL SRFTRNEFSL ESKSTIGVEF ATRSLNVNEK VIKAQIWDTA GQUERYRAITS AYYRGAVGAL LVDVTRHST FENVERWLRE LRDHTDPNIV VMLVGNKSDL RHLVAVQTED AKSFAENESL YFMETSALES TNVENAFSEV LTQIYHVVSK KAMEAGEDSG NVPSKGEKID VDVSAVKKTG CCSN
Source	Baculovirus
Target Names	RABA1D
Protein Names	Recommended name: Ras-related protein RABA1d Short name= AtRABA1d Alternative name(s): Ras-related protein Rab11B Short name= AtRab11B Ras-related protein SGBP Short name= AthSGBP
Expression Region	1-214
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.