



Recombinant Arabidopsis thaliana GDSL esterase/lipase 5 (GLIP5)

Product Code	CSB-BP874504DOA
Abbreviation	GLIP5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SSA7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	AKISH NGDNNVTALF LFGDSFLDAG NNNYINTTTL DQANFPPYGQ TFFGLPTGRF SDGRLISDFI AEYANLPLIP PFLEPGNSQK KLYGVNFASA GAGALVETFQ GSVINLRTQL DHYKKVERLW RTNFGKEESK KRISRAVYLI SIGSNDYSSI FLTNQSLPIS MSQHVDIVIG NLTTFIHEIY KIGGRKFGFL NVPDLGCFPA LRILQPKNDD SCLRASRLA SMHNRALTNL LFQMQRQVKG FKFSLFDMNK SLRLRMQHPS KFGFKEGEEA CCGTGKWRGV FSCGGKRIVK EYQLCENPKD YIFWDSLHLT QNTYNQFANL IWNGGHMSDS LVVGPYNINN LFQIP
Source	Baculovirus
Target Names	GLIP5
Protein Names	Recommended name: GDSL esterase/lipase 5 EC= 3.1.1.- Alternative name(s): Extracellular lipase 5
Expression Region	36-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.