



# Recombinant Tupaia belangeri Cellular tumor antigen p53 (TP53)

<b>Product Code</b>	CSB-BP874674TRC
<b>Abbreviation</b>	TP53
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9TTA1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Tupaia belangeri (Common tree shrew) (Tupaia glis belangeri)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEEPQSDPSV EPPLSQETFS DLWKLLPENN VLSPLPSQAM DDLMLSPDDI EQWFTEDPGP DEAPRMPEAA PPVAPAPAAP TPAAPAPAPS WPLSSSVPSQ KTYQGSYGFR LGFLHSGTAK SVTCTYSPDL NKLFCQLAKT CPVQLWVDSA PPPGRVRAM AIYKQSQYVT EVVRRCPHHE RCSDSDGLAP PQHLIRVEGN LHAEYSDDRN TFRHSVVVPY EPPEVGS DCT TIHYNMCNS SCMGGMNRRP ILTITLED SSGKLLGRNSF EVRICACPGR DRRTEENFR KKGESCPKLP TGSIKRALPT GSSSSPQPKK KPLDEEYFTL QIRGRERFEM LREINEALEL KDAMAGKESA GSRAHSSHLK SKKGQSTSRH RKL MFKTEGP DSD
<b>Source</b>	Baculovirus
<b>Target Names</b>	TP53
<b>Protein Names</b>	Recommended name: Cellular tumor antigen p53 Alternative name(s): Tumor suppressor p53
<b>Expression Region</b>	1-393
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.