



Recombinant Defective in cullin neddylation protein 1 (dcn-1)

Product Code	CSB-BP874714CXY
Abbreviation	dcn-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9U3C8
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MNRLKSDQKT KLRQFVQWTQ VTEAVSLNFL AKANWNIEYA MTLYFDNPNL FAGSTPQPSV DRSNIERLFN QYVDPKDKVG EKRMGPHGIN RLLTDLGYEA TDRRVLVLAW KFTAQTQCEF SLDEWVKGMT ALQADTVQNL RQRIDSINSG LESDKAKFHE LYLFAFNyak SAACRNLDLE TAICCWVDFL GQRSTIMTQW IDFLWAQENA AASRLAQNVG ASNAKQFKSV WISRDTWNLF WDFILLSKPD LSDYDDEGAW PVLIDQFVDY CRENLNPKP GNASNDQQME TPSYY
Source	Baculovirus
Target Names	dcn-1
Protein Names	Recommended name: Defective in cullin neddylation protein 1
Expression Region	1-295
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.