



Recombinant Human BTB/POZ domain-containing protein 10 (BTBD10)

Product Code	CSB-MP874798HU
Abbreviation	BTBD10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BSF8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAGRPHPYDG NSSDPENWDR KLHSRPRKLY KHSSTSSRIA KGGVDHTKMS LHGASGGHER SRDRRRSSDR SRDSSHERTE SQLTPCIRNV TSPTRQHHVE REKDHSSSRP SSPRPQKASP NGSISSAGNS SRNSSQSSSD GCKTAGEMV FVYENAKEGA RNIRTSERVV LIVDNTRFVV DPSIFTAQPN TMLGRMFGSG REHNFTRPNE KGEYEVAEGI GSTVFRAILD YYKTGIIRCP DGISIPELRE ACDYLCISFE YSTIKCRDLS ALMHELSNDG ARRQFEFYLE EMILPLMVAS AQSGERECHI VVLTDDDDVVD WDEEYPPQMG EEYSQIIYST KLYRFFKYIE NRDVAKSVLK ERGLKKIRLG IEGYPTYKEK VKKRPGRPE VIYNYVQRPF IRMSWEKEEG KSRHVDFQCV KSKSITNLAA AAADIPQDQL VVMHPTPQVD ELDILPIHPP SGNSDLDPDA QNPML
Source	Mammalian cell
Target Names	BTBD10
Protein Names	Recommended name: BTB/POZ domain-containing protein 10
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.