



# Recombinant Human BTB/POZ domain-containing protein 10 (BTBD10)

<b>Product Code</b>	CSB-YP874798HU
<b>Abbreviation</b>	BTBD10
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BSF8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAGRPHPYDG NSSDPENWDR KLHSRPRKLY KHSSTSSRIA KGGVDHTKMS LHGASGGHER SRDRRRSSDR SRDSSHERTE SQLTPCIRNV TSPTRQHHVE REKDHSSSRP SSPRPQKASP NGSISSAGNS SRNSSQSSSD GCKTAGEMV FVYENAKEGA RNIRTSERVV LIVDNTRFVV DPSIFTAQPN TMLGRMFGSG REHNFTRPNE KGEYEVAEGI GSTVFRAILD YYKTGIIRCP DGISIPELRE ACDYLCISFE YSTIKCRDLS ALMHELSNDG ARRQFEFYLE EMILPLMVAS AQSGERECHI VVLTDDDDVVD WDEEYPPQMG EEYSQIIYST KLYRFFKYIE NRDVAKSVLK ERGLKKIRLG IEGYPTYKEK VKKRPGRPE VIYNYVQRPF IRMSWEKEEG KSRHVDFQCV KSKSITNLAA AAADIPQDQL VVMHPTPQVD ELDILPIHPP SGNSDLDPDA QNPML
<b>Source</b>	Yeast
<b>Target Names</b>	BTBD10
<b>Protein Names</b>	Recommended name: BTB/POZ domain-containing protein 10
<b>Expression Region</b>	1-475
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.