



Recombinant Human TERF1-interacting nuclear factor 2 (TINF2)

Product Code	CSB-MP874799HU
Abbreviation	TINF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BSI4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ATPLVAGPA ALRFAAAASW QVVRGRCVEH FPRVLEFLRS LRAVAPGLVR YRHHERLCMG LKAKVVVELI LQGRPWAQVL KALNHHFPES GPIVRDPKAT KQDLRKILEA QETFYQQVKQ LSEAPVDLAS KLQELEQEYG EPFLAAMEKL LFEYLCQLEK ALPTPQAQQL QDVLSWMQPG VSITSSLAWR QYGVDMGWLL PECSVTDSVN LAEPMEQNPP QQQLRALHNP LPKAKPGTHL PQGPSSRTHP EPLAGRHFNL APLGRRRVQS QWASTRGGHK ERPTVMLFPF RNLGSPTQVI SKPESKEEHA IYTADLAMGT RAASTGKSKS PCQTLGGRAL KENPVDLPAT EQKENCLDCY MDPLRLSLLP PRARKPVCPP SLCSSVITIG DLVLDSDEEE NGQGEGKESL ENYQKTKFDT LIPTLCEYLP PSGHGAIPVS SCDCRDSSRP L
Source	Mammalian cell
Target Names	TINF2
Protein Names	Recommended name: TERF1-interacting nuclear factor 2 Alternative name(s): TRF1-interacting nuclear protein 2
Expression Region	2-451
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.