



Recombinant Human UPF0183 protein C16orf70 (C16orf70)

Product Code	CSB-EP874801HU-B
Abbreviation	C16orf70
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BSU1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLDLEVVPER SLGNEQWEFT LGMPLAQAVA ILQKHCRIIK NVQVLYSEQS PLSHDLILNL TQDGIKLMFD AFNQRLLKLVIE VCDLTKVKLK YCGVHFNSQA IAPTIEQIDQ SFGATHPGVY NSAEQLFHLLN FRGLSFSFQL DSWTEAPKYE PNFAHGLASL QIPHGATVVKR MYIYSGNSLQ DTKAPMMPLS CFLGNVYAES VDVLRDGTGP AGLRLRLLAA GCGPGLLADA KMRVFERSVY FGDSCQDVLS MLGSPHKVFY KSEDKMKIHS PSPHKQVPSK CNDYFFNYFT LGVDILFDAN THKVKKFVLH TNYPGHYNFN IYHRCEFKIP LAIKKENADG QTETCTTYSK WDNIQELLGH PVEKPVVLRH SSSPNNTNPF GSTFCFGLQR MIFEVMQNNH IASVTLYGPP RPGSHLRRTAE LP
Source	E.coli
Target Names	C16orf70
Protein Names	Recommended name: UPF0183 protein C16orf70
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.