



Recombinant Human Protein canopy homolog 3 (CNPY3)

Product Code	CSB-YP874802HU
Abbreviation	CNPY3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BT09
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GPSQAGAEEN DWVRLPSKCE VCKYVAVELK SAFEETGGTK EVIGTGYGIL DQKASGVKYT KSDLRLIEVT ETICKRLLDY SLHKERTGSN RFAKGMSETF ETLHNLVHKG VKVVMPIPYE LWNETSAEVA DLKKQCDVLV EEFEEVIEDW YRNHQEEDLT EFLCANHVLIK GKDTSCLAEQ WSGKKGDAA LGGKSKKKS SRAKAAGGRS SSSKQRKELG GLEGDPSPEE DEGIQKASPL THSPPEL
Source	Yeast
Target Names	CNPY3
Protein Names	Recommended name: Protein canopy homolog 3 Alternative name(s): CTG repeat protein 4a Expanded repeat-domain protein CAG/CTG 5 Protein associated with TLR4 Trinucleotide repeat-containing gene 5 protein
Expression Region	31-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.