



Recombinant Human THAP domain-containing protein 7 (THAP7)

Product Code	CSB-EP874803HU
Abbreviation	THAP7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BT49
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPRHCSAAGC CTRDTRETRN RGISFHRLPK KDNPRRGLWL ANCQRLDPSG QGLWDPASEY IYFCSKHFEEDCFELVIGISG YHRLKEGAVP TIFESFSKLR RTTKTKGHSY PPGPAEVSRL RRCRKRCSSEG RGPTTPFSPP PPADVTCFPV EEASAPATLP ASPAGRLEPG LSSPFSDLLG PLGAQADEAG CSAQPSPERQ PSPLEPRPVS PSAYMLRLPP PAGAYIQNEH SYQVGSALLW KRRAEAALDA LDKAQRQLQA CKRREQRLRL RLTKLQGERA REKRAQADAR QTLKEHVQDF AMQLSSSMA
Source	E.coli
Target Names	THAP7
Protein Names	Recommended name: THAP domain-containing protein 7
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.