



Recombinant Human REM2- and Rab-like small GTPase 1 (RSG1)

Product Code	CSB-MP874809HU
Abbreviation	RSG1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BU20
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MARPPVPGSV VVPNWHESAE GKEYLACILR KNRRRVFGLL ERPVLLPPVS IDTASYKIFV SGKSGVGKTA LVAKLAGLEV PVVHHETTGI QTTVVFWPAK LQASSRVVMF RFEFWDCGES ALKKFDHMLL AMENTDAFL FLFSFTDRAS FEDLPGQLAR IAGEAPGVVR MVIGSKFDQY MHTDVPERDL TAFRQAWELP LLRVKSVPGR RLADGRTLGD RAGLADVAHI LNGLAEQLWH QDQVAAGLLP NPPEAPE
Source	Mammalian cell
Target Names	RSG1
Protein Names	Recommended name: REM2- and Rab-like small GTPase 1
Expression Region	1-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.