



Recombinant Human Single-stranded DNA-binding protein 4 (SSBP4)

Product Code	CSB-MP874823HU
Abbreviation	SSBP4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BWG4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MYAKGGKGS A VPSDSQAREK LALYVY EYLL HIGAQKSAQT FLSEIRWEKN ITLGEPPGFL HSWWCVFWDL YCAAPDRREA CEHSGEAKAF QDYSAAAAPS PVMGSMAPGD TMAAGSMAAG FFQGPPGSQP SPHNP NAPMM GPHGQPFMSP RFPGGPRPTL RMPSQPPAGL PGSQPLLPGA MEPSRAQGH PSMGGPMQRV TPRGMASVG PQSYGGGMRP PPNSLAGPGL PAMNMGPVGR GPWASPSGNS IPYSSSSPGS YTGPPGGGGP PGTPI MPSPG DSTNSENMY TIMNPIGQGA GRANFPLGPG PEGPMAAMSA MEPHHVNGSL GSGDMDGLPK SSPGAVAGLS NAPGTPRDDG EMAAAGTFLH PFPSESYSPG MTMSV
Source	Mammalian cell
Target Names	SSBP4
Protein Names	Recommended name: Single-stranded DNA-binding protein 4
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.