



Recombinant Human Cell division cycle-associated protein 7 (CDCA7)

Product Code	CSB-MP874825HU
Abbreviation	CDCA7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BWT1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDARRVPQKD LRVKKNLKKF RYVKLISMET SSSSDDSCDS FASDNFANTR LQSVREGCRT RSQCRHSGPL RVAMKFPARS TRGATNKKAE SRQPSENSVT DSNSDSEDES GMNFLEKRAL NIKQNKAMLA KLMSELESFP GSFRGRHPLP GSDSQSRRPR RRTFPGVASR RNPERRARPL TRSRSRILGS LDALPMEEEE EEDKYMLVRK RKTVDGYMNE DDLPRSRRSR SSVTLPHIIR PVEEITEEEL ENVCSNSREK IYNRSLGSTC HQCRQKTIDT KTNCRNPDCW GVRGQFCGPC LRNRYGEEVR DALLDPNWHC PPCRGICNCS FCRQRDGRCA TGVLVYLAKY HGFGNVHAYL KSLKQEFEMQ A
Source	Mammalian cell
Target Names	CDCA7
Protein Names	Recommended name: Cell division cycle-associated protein 7 Alternative name(s): Protein JPO1
Expression Region	1-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene was identified as a c-Myc responsive gene, and behaves as a direct c-Myc target gene. Overexpression of this gene is found to enhance the transformation of lymphoblastoid cells, and it complements a transformation-defective Myc Box II mutant, suggesting its involvement in c-Myc-mediated cell transformation. Two alternatively spliced transcript variants encoding distinct isoforms have been reported.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.