



# Recombinant Human Oxysterol-binding protein-related protein 11 (OSBPL11)

<b>Product Code</b>	CSB-YP874829HU
<b>Abbreviation</b>	OSBPL11
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BXB4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>                     MQGGEPVSTM KVSESEGKLE GQATAVTPNK NSSCGGGISS SSSSRGGSAAK                      GWQYSDHMEN VYGYSLMKYTN LVTGWQYRFF VLNNEAGLLE                      YFVNEQSRNQ KPRGTLQLAG AVISPSDEDS HTFTVNAASG EQYKLRATDA                      KERQHWVSRL QICTQHHTEA IGKNNPPLKS RSFSLASSSN SPISQRRPSQ                      NAISFFNVGH SKLQSLSKRT NLPPDHLVEV REMMSHAEGQ QRDLIRRIEC                      LPTSGHLSSL DQDLLMLKAT SMATMNCND CFHILQLQHA SHQKGLPSG                      TTIEWLEPKI SLSNHYKNGA DQPFATDQSK PVAVPEEQPV AESGLLAREP                      EEINADDEIE DTCDHKEDDL GAVEEQRSVI LHLLSQLKLG MDLTRVVLPT                      FILEKRSLE MYADFMSPD LFIAITNGAT AEDRMIRFVE YYLTSFHEGR                      KGAIAKKPYN PIIGETFHCS WKMPKSEVAS SVFSSSSTQG VTNHAPLSGE                      SLTQVGSDCY TVRFVAEQVS HPPVSGFYA ECTERKMCVN AHVWTKSKFL                      GMSIGVTMVG EGILSLEHG EETYSLPCA YARSILTPW VELGGKSVN                      CAKTGYSASI TFHTKPFYGG KLHRVTAEVK HNITNTVVCV VQGEWNSVLE                      FTYSNGETKY VDLTKLAVTK KRVRPLEKQD PFESRRLWKN VTDSLRESEI                      DKATEHKHTL EERQRTEERH RTETGTPWKT KYFIKEGDGW VYHKPLWKII                      PTTQPAE                 </p>
<b>Source</b>	Yeast
<b>Target Names</b>	OSBPL11
<b>Protein Names</b>	Recommended name: Oxysterol-binding protein-related protein 11 Short name= ORP-11 Short name= OSBP-related protein 11
<b>Expression Region</b>	1-747
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a member of the oxysterol-binding protein (OSBP) family, a group of intracellular lipid receptors. Like most members, the encoded protein



contains an N-terminal pleckstrin homology domain and a highly conserved C-terminal OSBP-like sterol-binding domain.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.