



# Recombinant Human Inosine triphosphate pyrophosphatase (ITPA)

<b>Product Code</b>	CSB-YP874838HU
<b>Abbreviation</b>	ITPA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BY32
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AASLVGKKI VFVTGNAKKL EEVVQILGDK FPCTLVAQKI DLPEYQGEPD EISIQKCQEA VRQVQGPVLV EDTCLCFNAL GGLPGPYIKW FLEKLKPEGL HQLLAGFEDK SAYALCTFAL STGDPSQPVR LFRGRTSGRI VAPRGCQDFG WDPCFQPDGY EQTYAEMPKA EKNAVSHRFR ALLELQEYFG SLAA
<b>Source</b>	Yeast
<b>Target Names</b>	ITPA
<b>Protein Names</b>	Recommended name: Inosine triphosphate pyrophosphatase Short name= ITPase Short name= Inosine triphosphatase EC= 3.6.1.19 Alternative name(s): Non-canonical purine NTP pyrophosphatase Non-standard purine NTP pyrophosphatas
<b>Expression Region</b>	2-194
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein hydrolyzes inosine triphosphate and deoxyinosine triphosphate to the monophosphate nucleotide and diphosphate. The encoded protein, which is a member of the HAM1 NTPase protein family, is found in the cytoplasm and acts as a homodimer. Defects in the encoded protein can result in inosine triphosphate pyrophosphorylase deficiency. Two transcript variants encoding two different isoforms have been found for this gene. Also, at least two other transcript variants have been identified which are probably regulatory rather than protein-coding.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.