



Recombinant Human UBX domain-containing protein 6 (UBXN6)

| | |
|--------------------------|--|
| Product Code | CSB-MP874859HU |
| Abbreviation | UBXN6 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9BZV1 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MKKFFQEFKA DIKFKSAGPG QKLKESVGEK AHKEKPNQPA PRPPRQGPTN EAQMAAAAAL ARLEQKQSRA WGPTSQDTIR NQVRKELQAE ATVSGSPEAP GTNVVSEPRE EGSAHLAVPG VYFTCPLTGA TLRKDQRDAC IKEALLHFS TDPVAASIMK IYTFNKDQDR VKLGVDTIAK YLDNIHLHPE EEKYRKIKLQ NKVLFQERINC LEGTHEFFEA IGFQKVLLPA QDQEDPEEFY VLSETTLAQP QSLERHKEQL LAAEPVRAKL DRQRRVFQPS PLASQFELPG DFFNLTAEEI KREQRLRSEA VERLSVLRTK AMREKEEQRG LRKYNYTLR VRLPDGCLLQ GTFYARERLG AVYGFVREAL QSDWLPFELL ASGGQKLSER ENLALNECGL VPSALLTFSW DMAVLEDIKA AGAEPDSILK PELLSAIEKL L |
| Source | Mammalian cell |
| Target Names | UBXN6 |
| Protein Names | Recommended name: UBX domain-containing protein 6 Alternative name(s): UBX domain-containing protein 1 |
| Expression Region | 1-441 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.