



Recombinant Human Forkhead box protein Q1 (FOXQ1)

Product Code	CSB-EP874860HU
Abbreviation	FOXQ1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C009
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKLEVFVPRA AHGDKQGSDL EGAGGSDAPS PLSAAGDDSL GSDGDCAANS PAAGGGARDT QGDGEQSAGG GPGAEEAIPA AAAAVVAEG AEAGAAGPGA GGAGSGEGAR SKPYTRRPKP PYSYALIAM AIRDSAGGRL TLAEINEYLM GKFPFFRGSY TGWRNSVRHN LSLNDCFVKV LRDPSRPWGK DNYWMLNPNS EYTFADGVFR RRRKRLSHRA PVPAPGLRPE EAPGLPAAPP PAPAAPASPR MRSPARQEER ASPAGKFSSS FAIDSILRKP FRSRRLRDTA PGTTLQWGAA PCPPLPAFPA LLPAAPCRAL LPLCAYGAGE PARLGAREAE VPPTAPLLL APLPAAAPAK PLRGPAAGGA HLYCPLRLPA ALQAASVRRP GPHLPYPVET LLA
Source	E.coli
Target Names	FOXQ1
Protein Names	Recommended name: Forkhead box protein Q1 Alternative name(s): HNF-3/forkhead-like protein 1 Short name= HFH-1 Hepatocyte nuclear factor 3 forkhead homolog 1
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.