



Recombinant *Schizosaccharomyces pombe* RuvB-like helicase 1 (rvb1)

Product Code	CSB-MP874872SXV
Abbreviation	rvb1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C0X6
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MVQISEVKGN GRDNRITTHS HIKGLGLKED GTCESVGGGF IGQEKAREAC GIITDLIKSK KFGGKGVLF A GGAGTGKTAL ALAIAQELGP KVPFCPMVGS EVYSSEIKKT EALMENFRRA IGLRVKETKE VYEGETEMV PEEAENPLGG YGKTISHVLL GLKTHKGTKQ LKLDPSIYES LQREQVSTGD VIYIEANTGA VKRVGRSDAY ATEFDLEAEE YVPMKGEVH KRKEIVQDVT LHDLDIANAR PQGGQDIMSM MGQLMKPKKT EITDKLRGEI NKVVNKYIEQ GIAELIPGVL FIDEVHMLDI ECFTYLNQAL ESTISPIVIF ASNRGICTIR GTEDIQAPHG IPTDLLDRLL IVRTLPHYSES EIRSILQIRA KVENIILTDE CLDKLAQEGS RTSLRYVIQL LTPVSIASL HGNKEIGVQD IEECNDLFLD ARRSAQVVKS SSGFLQ
Source	Mammalian cell
Target Names	rvb1
Protein Names	Recommended name: RuvB-like helicase 1 EC= 3.6.4.12
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.