



Recombinant Schizosaccharomyces pombe

Uncharacterized protein PB24D3.07c

(SPAPB24D3.07c)

Product Code	CSB-EP874873SXV-B
Abbreviation	SPAPB24D3.07c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C0Y7
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	ANASNSSV PTPDNTLVVS YTNTSCYTSG PLSPDTRFNR TTRGTFSKVR DALKFRLNGP IHHWDIANEL FDTALGVEII DTQYGINNRT SRDWCTAVSA LEKGDLEIFA AAFTAFDDVN PTKEVVPDAL VGTLALWAKY SYEYDLTSIV SLFGKNWTV D ELGWGTYIAH SLNEISSNTT GAANATLFID TQTKSDCYKL ASTIESWKYA PLSAFPNYGP FYIYGQCVAT FTSGYSPLVE PAFTFASALN TTLNSFPNGT LPTQTQVVGD IGVKFLQYFA
Source	E.coli
Target Names	SPAPB24D3.07c
Protein Names	Recommended name: Uncharacterized protein PB24D3.07c
Expression Region	23-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.