



Recombinant Arabidopsis thaliana Agamous-like MADS-box protein AGL86 (AGL86)

Product Code	CSB-EP874905DOA-B
Abbreviation	AGL86
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C6V3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MRSKIKLSLI ANKTSRRTTF RKRKGGITNK LHELTTLCGV KACAVISSPY ENPVVWPSTE GVQEAVSMFM ERPATEQSKL MMSHETYLQD KITKETKKLE SLRRENRESQ LRQFMFDCVE GKMSEHQYGA RDLQDLSLYI DHYINQLNSS VMLLTNNGAS SSSFPPPLHT SVAGAGAGAG AAPLVVAGAG AAPLAVAGAG ASPLAVAGVG AAPLAVAGAG PPMAQNQYEP IQPYIPTAFS DNIQYQAPVD FNHQIQHGIY DNLSLDPNHQ YPFQDDPFME MLMEYPYEQV GYAAEHAHIP FMNGNYNYH QPPTVGLTTT GHMPSNNATT TTTTNTTV
Source	E.coli
Target Names	AGL86
Protein Names	Recommended name: Agamous-like MADS-box protein AGL86
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.