



Recombinant Arabidopsis thaliana Putative F-box protein At3g10430 (At3g10430)

Product Code	CSB-BP874937DOA
Abbreviation	At3g10430
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9CAE7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGSSLPFDLI LEILQRTPAE SLLRFKSTCK KWYELISNDK RFMYKHLDKS TKRFLRIENR ERVQILD PVT EILAVSTIPN ELRHKYFTLI HCDGLMLGMC YEELGSDPNL AVWNPV MRKI KWIKSPPLV CYWGSDYLG F GYDKTFRDNY KILRFTYLGD DDDDESYPKC QIYEFNSGSW RSIEAKFDGE IDVEVDGVS V NGSMYWIELQ EKKNFILSFD FSKETFNRIC DSPLYWDIKR LGCFNGDRLS LLQQNEATRE IHVSVTNKLS DEVVSFSRYF SVTLPDIPLF RSIDILPLPG CFIGKNRNII AWCEQLVGKY DDHTSII LYE IGEGGVTKRI RTGRHEADED YNDLDNCSYV YVPSLIPVPE
Source	Baculovirus
Target Names	At3g10430
Protein Names	Recommended name: Putative F-box protein At3g10430
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.