



# Recombinant Arabidopsis thaliana Probable WRKY transcription factor 36 (WRKY36)

<b>Product Code</b>	CSB-BP874941DOA
<b>Abbreviation</b>	WRKY36
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9CAR4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MIKEETVSYF QTFDGVMAES DKEEELDATK AKVEKVREEN EKLKLLLSTI LNNYNLSLQMQ VSKVLGQQQG ASSMELDHID RQDENNDYDV DISLRLGRSE QKISKKEENK VDKISTKNVE ESKDKRSALG FGFQIQSYEA SKLDDLRCRQV KLANAENKCV SSRKDVKSVR NENHQDVLEE HEQTGLKKTR VCVKASCEDP SINDGCQWRK YGQKTAKTNP LPRAYYRCSM SSNCPVRKQV QRCGEEETSA FMTTYEGNHD HPLPMEASHM AAGTSA AASL LQSGSSSSSS STSASLSYFF PFHHFSISTT NSHPTVTLDL TRPNYPNQLP DDYPLSSSSF SLNFSSPDPP PPSSHDHTLN FSGLRTQAPL STDSLLARYR TRLSGQQ
<b>Source</b>	Baculovirus
<b>Target Names</b>	WRKY36
<b>Protein Names</b>	Recommended name: Probable WRKY transcription factor 36 Alternative name(s): WRKY DNA-binding protein 36
<b>Expression Region</b>	1-387
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.