



# Recombinant Mouse ATP synthase subunit b, mitochondrial (Atp5f1)

<b>Product Code</b>	CSB-YP875117MO
<b>Abbreviation</b>	Atp5f1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9CQQ7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	PLPPLPEY GGKVRGLIP EEFFQFLYPK TGVTGPYVLG TGLSLYFLSK EIVVITPETF STISVVGLIV YVIKKYGASF GEFIDKLNEE KIAQLEEVKQ SSMKQIQDAI DMEKAQQALV QKRHYLFDVQ RNNIALALEV TYRERLHKAY KEVKNRLDYH ISVQNMRRK EEEHMIDWVE KHVVKISISVQ QEKETIAKCI EDLKLLAKKA QAQPIM
<b>Source</b>	Yeast
<b>Target Names</b>	Atp5f1
<b>Protein Names</b>	Recommended name: ATP synthase subunit b, mitochondrial Short name= ATPase subunit b
<b>Expression Region</b>	43-256
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a subunit of mitochondrial ATP synthase. Mitochondrial ATP synthase catalyzes ATP synthesis, utilizing an electrochemical gradient of protons across the inner membrane during oxidative phosphorylation. ATP synthase is composed of two linked multi-subunit complexes: the soluble catalytic core, F1, and the membrane-spanning component, Fo, comprising the proton channel. The catalytic portion of mitochondrial ATP synthase consists of 5 different subunits (alpha, beta, gamma, delta, and epsilon) assembled with a stoichiometry of 3 alpha, 3 beta, and a single representative of the other 3. The proton channel seems to have nine subunits (a, b, c, d, e, f, g, F6 and 8). This gene encodes the b subunit of the proton channel.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.