



Recombinant Mouse WD repeat, SAM and U-box domain-containing protein 1 (Wdsub1)

Product Code	CSB-BP875173MO
Abbreviation	Wdsub1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9D0I6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MVRLIHTLAD HGDDVSCCAF SAALLATCSL DKTIRLYSLS DFAELPYSP KFHTYAVHCC CFSPSGHVLA SCSTDGTTVL WSSHSGHTLT VLEQPGGSPV RVCCFSPDSA YLASGAADGS IALWNAQTYK LYRCGSVKDS SLVACAFSPD GGLFVTGSSG GDLTVWDDRM RCLHSEKAHD LGITCCSFSS QPLSGGEGLO SYQLASCGQD CEIKLWAVTI TRVLGFELKY KSTLSGHCAP VLACAFSHDG KMLASGSVDK SVIIHGIGPQ SVLHTLTQHT RYVTTCAFAP NTLLLATGSM DKTVNIWQFD LETPCQGS MN DPLKHFTTEW SEEDVSVWLR AQGLEDLVGI FRANNIDGKE LLHLTKESLA GDLKIESLGL RSKVLRISIE LRAKMDSLSS GIPDEFICPI TRELMKDPVI ASDGYSYERE AMESWIHKKK RTSPMTNLAL PSLVLT PNRT LKMAINRWLE THEK
Source	Baculovirus
Target Names	Wdsub1
Protein Names	Recommended name: WD repeat, SAM and U-box domain-containing protein 1
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.