



# Recombinant Mouse 28S ribosomal protein S25, mitochondrial (Mrps25)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP875178MO   |
| <b>Abbreviation</b>      | Mrps25   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q9D125   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Mus musculus (Mouse)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MPMKGRFPIR RTLQYLGRGD VVKESVKIM TVNYNTYGEL GEGARKFVFF<br>NIPQIQYKNP WVQIMMFKNM TPSPFLRFYL DSGEQVLVDV ETKSNKEIME<br>HIKKILGKKE ETLREEELEK QQRFHGPNFG PRKYCLRECM CEVEGQVPCP<br>GLVPLPKEMT GKYKAALKAS T   |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | Mrps25   |
| <b>Protein Names</b>     | Recommended name: 28S ribosomal protein S25, mitochondrial Short name=MRP-S25 Short name= S25mt  |
| <b>Expression Region</b> | 1-171  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Target Details</b>    | Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein. A pseudogene corresponding to this gene is found on chromosome 4. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final   |



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.