



Recombinant Mouse Neuroguidin (Ngdn)

Product Code	CSB-YP875269MO
Abbreviation	Ngdn
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9DB96
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	AAPEVLESD VSSITLLKN LQEQVMAVTA QIQALTTKVR AGTYSTEKGL SFLEVKDQLL LMYLMDLSHL ILDKASGASL QGHPAVLRLV EIRTVLEKLR PLDQKLKYQI DKLVKTAVTG SLENDPLRF KPHPSNMVSK LSSEDEESEE AEEGQSEASG KKSAGSACK YVPPRLVPVH YDETEAEREQ KRLEKAKRRA LSSSVIRELK EQYSDAPEEI RDARHPHVTR QSQEDQHRVN YEESMMVRLS VSKREKGLRR RASAMSSQLH SLTHFSDISA LTGGTAHLDE DQNPVKKRKK LPKKGRKKKG FRRRW
Source	Yeast
Target Names	Ngdn
Protein Names	Recommended name: Neuroguidin Alternative name(s): EIF4E-binding protein
Expression Region	2-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>