



Recombinant Mouse Alpha-parvin (Parva)

Product Code	CSB-MP875329MO
Abbreviation	Parva
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9EPC1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	ATSPQKSPL VPKSPTPKSP PSRKKDDSFL GKLGGTLARR KKAKEVSEFQ EEGMNAINLP LSPISFELDP EDTLLEENEV RTMVDPNRN DPKLQELMKV LIDWINDVLV GERIIVKDLA EDLYDGQVLQ KLFEKLESEK LNVAEVTQSE IAQKQKLQTV LEKINETLKL PPRSİKWNVD SVHAKNLVAI LHLLVALSQY FRAPIRLPDH VSIQVVVVQK REGILQSRQI QEEITGNTA LSGRHERDAF DTLFDHAPDK LNVVKKTLIT FVNKHLNKLN LEVTELETQF ADGVYLVLLM GLLEGYFVPL HSFFLTPDSF EQKVLNVSFA FELMQDGGLE KPKRPEDIV NCDLKSTLRV LYNLFTKYRN VE
Source	Mammalian cell
Target Names	Parva
Protein Names	Recommended name: Alpha-parvin Alternative name(s): Actopaxin
Expression Region	2-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>