



Recombinant Mouse Probable folate receptor delta (Folr4)

Product Code	CSB-YP875340MO
Abbreviation	Folr4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9EQF4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	G DKLLSVCMNS KRHKQEPGPE DELYQECRPW EDNACCTRST SWEAHLEEPL LFNFSMMHCG LLTPACRKHF IQAICFHECS PNLGPWIQPV VPNGQEEQRV WGVPLCQEDC EDWWRACHSS LTCKSNWLHG WDWSEEKHC PAHEPCLPFS YHFPTDDLK EKIWNNTFKA SPERRNSGRC LQKWFPTLS NPNVEVALHF AGSALAPQLS YTLPAFSLCL LFHP
Source	Yeast
Target Names	Izumo1r
Protein Names	Recommended name: Probable folate receptor delta Short name= FR-delta Alternative name(s): Folate receptor 4 Folate-binding protein 3
Expression Region	20-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.