



Recombinant Arabidopsis thaliana Pumilio homolog 26 (APUM26)

Product Code	CSB-BP875439DOA
Abbreviation	APUM26
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FGE7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MLNLDFTHKT TQATPRLHAV ATEFLRVSND VAEHLKLSSK LTSDPYLFVE FVKTIRGFLS VQTALGLSGE IDTVFLQVIK GWFPDLITET FSFLIVVRII NLFNKRANSK VYPDILRRIG NNALYLTRNP LRGICLVEKA INVRDPDCTV FIALKLHSHY VELSFEEELGS NIVEKLLSVG ESGICGV
Source	Baculovirus
Target Names	APUM26
Protein Names	Recommended name: Pumilio homolog 26 Short name= APUM-26 Short name= AtPUM26
Expression Region	1-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.