



Recombinant Arabidopsis thaliana F-box protein At5g51370 (At5g51370)

Product Code	CSB-MP875442DOA
Abbreviation	At5g51370
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FGN4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSEFDKHVQN HIKNEKPFNS IKPSEIDSTL SLSDSLLLKI IEKLPESQSN DVSLVCKRWL NLQGQRLRSL KLLDFDFLLS ERLTTRFPNL THVDLVNACM NPRVNSGILF CHKSISFHLS SDSSNWEFLE ENLLHSDVID RGLRILSRES FDLLNLKVIN ASELGLLSLA GDCSDLQELE LHKCNDNLLH GIAACKNLKG LRLVGSVDGL YSSSVSDIGL TFLAQGCRSL VKLELSGCEG SFDGIKAIGQ CCEVLEELSI CDHRMDDGWI AALSYFESLK ILRISSCRKI DASPGPEKLL RSCPAMESLQ LKRCCLNDKE GIKALFKVCD GATEVNIQDC WGLSDDCFSL AKAFRRVRFL SLEGCSVLTS GGLESVILHW EELESMRVVS CKSIKDSEIS PALSSLFSLI KELTWRPDTR SHLSSSLEGA GIGIRGSKFF KKSVLH
Source	Mammalian cell
Target Names	At5g51370
Protein Names	Recommended name: F-box protein At5g51370
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.