



Recombinant Arabidopsis thaliana B3 domain-containing protein At5g60140 (At5g60140)

Product Code	CSB-YP875449DOA
Abbreviation	At5g60140
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FHB0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNTRGNYSNG FTSKFFKPYL PSESGDDLVL PISFNSCLPK SLPKTVTVRS ISGNIWKLGF KKCGGEVERF VMVSGWKKIV RDENLNGGDL LSFEDGSHF FNFSIFDHET TCKRLKRSSE QSKDIIKVG S DCEEESQASD DVIVLNSDDS DDSDNDYSVE DDNVAEDDDG LEDEV DVEAE DGYDAKDSDG LEDEDDDEAE DGYDAKDDDG LEDEDDLEDE DDERRYLDDH ENPYFTMTLN PKKKSQ LHIP AHVIRDYDLT FPERITVVDE MGALEKAIKI QKNGCIFVKG FGSFIRRNKM KMTDKMICEL KRTGGNLVHT IKVKIISG
Source	Yeast
Target Names	At5g60140
Protein Names	Recommended name: B3 domain-containing protein At5g60140
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.