



Recombinant Arabidopsis thaliana Calcium-binding protein CML37 (CML37)

Product Code	CSB-BP875462DOA
Abbreviation	CML37
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FIH9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MTLAKNQKSS LSRLYKKVSS KRSESSRNLE DESRTSSNSS GSSSLNVNEL RTVFDYMDAN SDGKISGEEL QSCVSLGGGA LSSREVEEVV KTSVDVDGDGF IDFEEFLKLM EGEDGSDEER RKELKEAFGM YVMEGEEFIT AASLRRTLSR LGESCTVDAC KVMIRGFDQN DDGVLSEDEF VLMMR
Source	Baculovirus
Target Names	CML37
Protein Names	Recommended name: Calcium-binding protein CML37 Alternative name(s): Calmodulin-like protein 37
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.