



Recombinant Human Guanine nucleotide-binding protein G (olf) subunit alpha (GNAL)

Product Code	CSB-MP009592HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P38405
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGCLGGNSKT TEDQGVDEKE RREANKKIEK QLQKERLAYK ATHRLLLLGA GESGKSTIVK QMRILHVNGF NPEEKKQKIL DIRKNVKDAI VTIVSAMSTI IPPVPLANPE NQFRSDYIKS IAPITDFEYS QEFFDHVKKL WDDEGVKACF ERSNEYQLID CAQYFLERID SVSLVDYTPPT DQDLLRCRVL TSGIFETRFQ VDKVNFMFD VGGQRDERRK WIQCFNDVTA IYVAACSSY NMVIREDNNT NRLRESLDF ESIWNNRWLR TISIILFLNK QDMLAEKVLA GKSKIEDYFP EYANYTVPED ATPDAGEDPK VTRAKFFIRD LFLRISTATG DGKHYCYPHF TCAVDTENIR RVFNDCRDII QRMHLKQYEL L
Source	Mammalian cell
Target Names	GNAL
Protein Names	Recommended name: Guanine nucleotide-binding protein G(olf) subunit alpha Alternative name(s): Adenylate cyclase-stimulating G alpha protein, olfactory type
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.