



Recombinant Arabidopsis thaliana Agamous-like MADS-box protein AGL62 (AGL62)

Product Code	CSB-EP875481DOA
Abbreviation	AGL62
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FKK2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MVKKSKGRQK IEMVKMKNES NLQVTFKRR SGLFKKASEL CTLCGAEVAI VVFPGRKVF SFGHPNVDSV IDRFINNPL PPHQHNNMQL RETRRNSIVQ DLNNHLTQVL SQLETEKKKY DELKKIREKT KALGNWWEDP VEELALSQLE GFKGNLENLK KVVTVESRF FQANVPNFYV GSSSNNAAFG IDDGSHINPD MDLFSQRRMM DINAFNYNQN QIHPNHALPP FGNNAYGINE GFVPEYNVNF RPEYNPNQIQ IQNQNVQIQ IQNQSFKREN ISEYEHGGY PPQSRSDYY
Source	E.coli
Target Names	AGL62
Protein Names	Recommended name: Agamous-like MADS-box protein AGL62
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.