



Recombinant Arabidopsis thaliana Zinc finger CCCH domain-containing protein 59 (ZFWD3)

Product Code	CSB-EP875483DOA-B
Abbreviation	ZFWD3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FKR9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGHQSSWMKQ TVIDSVVQSS KVGVFMSMFSM DYKAPRKYSH GRNFGVARQQ DFAADIVTRR PYAPYDNGMK KGPNKLSRNL VWTSKEYKSP EGNRPRQNAA NGSAPQVIG TGHRVSNQPR KNAVYGRSS SLSDTRGCGP RLNGSPKKS SV CNFWKDGNC KGEKCQFLHS WSCFPGLAMV AALEGHKNDI KGIALPQGS D KLFSVSGDGT LLIWDCNSGQ CVRSINLQAE AGSLISEGPW VFLGLPNAVK AFNVQNSKDV HLEGVVGQVH AMTAANGMLF AGTSSGSILV WKATDSESDP FKYLTSLEGH HSGEVTCFVV GGEVLYSGSV DKTIKVWDLN TLQCRMTLKQ HIGTVTSLLC WDKCLISSSL DGTIKLWACS ENESLKVVQT RKQELSVHTL CGMHDAEAKP IMFCSYQNGA VGIFDLPSFE ERGKMFSTQT ICTLTIGPGG LLSFGDKSGN LRVWVSLASGT KV
Source	E.coli
Target Names	ZFWD3
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 59 Short name= AtC3H59 Alternative name(s): Zinc finger CCCH domain and WD40 repeat-containing protein 3
Expression Region	1-472
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.