



Recombinant Arabidopsis thaliana Putative fatty acyl-CoA reductase 7 (FAR7)

Product Code	CSB-BP875508DOA
Abbreviation	FAR7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FMQ9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEPNCVQFLE NRTILVTGAS GFLAKVLVER ILRLQPNVKR LYLLVRASDK KSAEQRYDVA LGINTFGAIN VLNFAKKCVK PKLLLVSTV YVCGERPGHI VEKHFAMGES LNGKNKVDIN TERRLADQKS KQFKEQGCSE EETEAMKDF GLKRARLYGW PNTYVFTKAM GEMLLGHYRE TMPVIRPT IITSTFSDPF PGWIEGLKTV DSVIFYGKG ILKCFVLDQK TVCDIIPVDM VVNAMIAIAA DHCHDSGSHT VYHVGSSNQN PVIYKQIYEM MSRYFMKSPL VGRNGMLIVP KVTRISTLAR FRVYTNLRYK LPIQILGLLS VISLSQRDKF ALHNRKFKMA MRLVKLYKPY VLFKGFDDK NMETLRIKNE AKDMEKLFGT NPKCIDWEDY FMIRISLAS
Source	Baculovirus
Target Names	FAR7
Protein Names	Recommended name: Putative fatty acyl-CoA reductase 7 EC= 1.2.1.n2
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.