



Recombinant Arabidopsis thaliana Flavin-containing monooxygenase FMO GS-OX-like 2 (At1g12200)

Product Code	CSB-EP875548DOA
Abbreviation	At1g12200
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FWW9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MATSHPDPTT SRHVAVIGAG AAGLVAAREL RREGHSVVVL ERGSQIGGVW AYTSQVEPDP LSLDPTRPVV HSSLYRSLRT NIPRECMGFT DFPFATRPHD GSRDPRRHPA HTEVLAYLRD FAKEFDIEEM VRFETEVVKA EQVAAEGEER GKWRVESRSS DGVVDEIYDA VVVCNGHYTE PRHALITGID SWPGKQIHSH NYRVPDQFKD QVVIVIGSSA SGVDICRDIA QVAKEVHVSS RSTSPDYEK LTGYENLWLH STIQIAREDG SVVFENGKTI YADTIMHCTG YKYYFPFLDT KGEVTVDDNR VGPLYKHVFP PALAPSLFI GLPWQITPPF MFELQSKWVA AVLSGRVSLP SQDEMMDTK AFYDKLEASG IPKRYTHLMP DDSQFEYDNW LADQCEYPRI EKWREQMFYI GFKRIYAQSS TYRDNWDDDH LIVEAYDDFV KFMSSYQELL PMLKT
Source	E.coli
Target Names	At1g12200
Protein Names	Recommended name: Flavin-containing monooxygenase FMO GS-OX-like 2 EC= 1.8.-.- Alternative name(s): Flavin-monooxygenase glucosinolate S- oxygenase-like 2
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.