



Recombinant Arabidopsis thaliana Vacuolar protein sorting-associated protein 20 homolog 2 (VPS20.2)

Product Code	CSB-EP875557DOA-B
Abbreviation	VPS20.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FY89
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	GNLFVKKPK ITEVDRAILS LKTQRRKLGQ YQQQLEKVIE AEKQAARDLI REKRKDRALL ALKKKRTQEE LLKQVDQWLI NVEQQLADIE LTSKQKAVFE SLKQGNNAIK AIQSEVNLDL VQKLMDDTAE AKAYQDELSA ILGEKLSAED EEEILAEFDN LESLLIVEDM PEVPTTELMP EEPEKMDLPD VPTKAPVASN ETTSTKRKVL EEPLA
Source	E.coli
Target Names	VPS20.2
Protein Names	Recommended name: Vacuolar protein sorting-associated protein 20 homolog 2 Short name= AtVPS20-2 Alternative name(s): Charged multivesicular body protein 6 homolog 2 ESCRT-III complex subunit VPS20 homolog 2
Expression Region	2-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.