



# Recombinant Arabidopsis thaliana Putative pentatricopeptide repeat-containing protein At1g02420 (At1g02420)

<b>Product Code</b>	CSB-MP875562DOA
<b>Abbreviation</b>	At1g02420
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9FZ19
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MMMILKPLSS HHVSNFRLSV SFLHSVALSD AKVPVEEEGD DAETVFRMIN GSNLQVELKE SLSSSGIHLS KDLIDRVLKR VRFSHGNPIQ TLEFYRYASA IRGFYHSSFS LDTMLYILGR NRKFDQIWEL LIETKRKDRS LISPRTMQVV LGRVAKLCSV RQTVESFWKF KRLVPDFFDT ACFNALLRTL CQEKSMTDAR NVYHSLKHQF QPDLQTFNIL LSGWKSSEEA EAFFEEMK GK GLKPDVVVTYN SLIDVYCKDR EIEKAYKLID KMREEEETPD VITYTTVIGG LGLIGQPDKA REVLKEMKEY GCYPDVAAYN AAIRNFCIAR RLGADADKLVD EMVKKGLSPN ATTYNLFFRV LSLANDLGRS WELYVRMLGN ECLPNTQSCM FLIKMFKRHE KVDAMAMLWE DMVVKGFGSY SLVSDVLLDL LCDLAKVEEA EKCLLEMVEK GHRPSNVSEFK RIKLLMELAN KHDEVNLIQ KMAIFSTEIP R
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At1g02420
<b>Protein Names</b>	Recommended name: Putative pentatricopeptide repeat-containing protein At1g02420
<b>Expression Region</b>	1-491
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.