



# Recombinant *Macaca fascicularis* Calpain-1 catalytic subunit (CAPN1), partial

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP875604MOV   |
| <b>Abbreviation</b>      | CAPN1   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9GLG2  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Macaca fascicularis</i> (Crab-eating macaque) ( <i>Cynomolgus</i> monkey)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | CAPN1   |
| <b>Protein Names</b>     | Recommended name: Calpain-1 catalytic subunit EC= 3.4.22.52 Alternative name(s): Calcium-activated neutral proteinase 1 Short name= CANP 1 Calpain mu-type Calpain-1 large subunit Micromolar-calpain Short name= muCA  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Partial   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |