



Recombinant Human Tubulointerstitial nephritis antigen-like (TINAGL1)

Product Code	CSB-MP875630HU
Abbreviation	TINAGL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZM7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AQQGRGRRE LAPGLHLRGI RDAGGRYCQE QDLCCRGRAD DCALPYLGAI CYCDLFCNRT VSDCCPDFWD FCLGVPPPPF PIQGCMHGGR IYPVLGTYWD NCNRCTCQEN RQWQCDQEPC LVDPDMIKAI NQGNYGWQAG NHSAFWGRTL DEGIRYRLGT IRPSSSVMMN HEIYTVLNPG EVLPTAFEAS EKWPNIHEP LDQGNCAWS AFSTAAVASD RVSISLGHM TPVLSPQNL SCDTHQQQGC RGGRLDGAWW FLRRRGVVSD HCYPFSGRER DEAGPAPPCM MHSRAMGRGK RQATAHCPNS YVNNNDIYQV TPVYRLGSND KEIMKELMEN GPVQALMEVH EDFFLYKGGI YSHTPVSLGR PERYRRHGTH SVKITGWGEE TLPDGRTLKY WTAANSWGPA WGERGHFRIV RGVNECDIES FVLGVWGRVG MEDMGHH
Source	Mammalian cell
Target Names	TINAGL1
Protein Names	Recommended name: Tubulointerstitial nephritis antigen-like Alternative name(s): Glucocorticoid-inducible protein 5 Oxidized LDL-responsive gene 2 protein Short name= OLRG-2 Tubulointerstitial nephritis antigen-related protein S
Expression Region	22-467
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.