



# Recombinant Human Poly [ADP-ribose] polymerase 12 (PARP12)

<b>Product Code</b>	CSB-BP875642HU
<b>Abbreviation</b>	PARP12
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H0J9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAQAGVVGEV TQVLCAAGGA LELPELRRRL RMGLSADALE RLLRQRGRFV VAVRAGGAAA APERVVLAAS PLRLCRAHQG SKPGCVGLCA QLHLCRFMVY GACKFLRAGK NCRNSHSLTT EHNLSVLRTH GVDHLSYNEL CQLLFQNDPW LLPEICQHYN KGDGPHGSCA FQKQCIKLI CQYFLQGECK FGTSCKRSHD FSNSENLEKL EKLGMSSDLV SRLPTIYRNA HDIKNKSSAP SRVPPLFVPQ GTSERKDSSG SVSPNTLSQE EGDQICLYHI RKSCSFQDKC HRVHFHLPYR WQFLDRGKWE DLDNMELIEE AYCNPKIERI LCSESASTFH SHCLNFNAMT YGATQARRLS TASSVTKPPH FILTDDWIWY WSDEFGSWQE YGRQGTVHPV TTVSSSDVEK AYLAYCTPGS DGQAATLKFQ AGKHNYELDF KAFVQKNLVY GTTKKVCRRP KYVSPQDVTT MQTCNTKFPG PKSIPDYWDS SALPDPGFQK ITLSSSSEY QKVWNLFNRT LPFYFVQKIE RVQNLALWEV YQWQKGQMOK QNGGKAVDER QLFHGTAIF VDAICQQNFD WRVCGVHGTS YGKGSYFARD AAYSHHYSKS DTQHTMFLA RVLVGEFVRG NASFVRPPAK EGWSNAFYDS CVNSVSDPSI FVIFEKHQVY PEYVIQYTTS SKPSVTPSIL LALGSLFSSR Q
<b>Source</b>	Baculovirus
<b>Target Names</b>	PARP12
<b>Protein Names</b>	Recommended name: Poly [ADP-ribose] polymerase 12 Short name= PARP-12 EC= 2.4.2.30 Alternative name(s): Zinc finger CCCH domain-containing protein 1
<b>Expression Region</b>	1-701
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.