



Recombinant Human Polyadenylate-binding protein 3 (PABPC3)

| | |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-EP875663HU |
| Abbreviation | PABPC3 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9H361 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MNPSTPSYPT ASLYVGD LHP DVTEAMLYEK FSPAGPILSI RICRDLITSG SSNYAYVNFQ HTKDAEHALD TMNFDVIK GK PVRIMWSQRD PSLRKSGVGN IFVKNLDKSI NNKALYDTVS AFGNILSCNV VCDENGSKGY GFVHFETHEA AERAIKKMNG MLLNGRKFV GQFKSRKERE AELGARAKEF PNVYIKNFGE DMDDERLKDL FGKFGPALS V KVM TDESGKS KGF GFVS FER HEDA QKAVDE MNGKELNGKQ IYVGRAQKV ERQTELKRTF EQMKQDRITR YQVVNLYVKN LDDGIDDERL RKA FSPFGTI TSAKVMMEGG RSKGFGFVCF SSPEEATKAV TEMNGRIVAT KPLYVALAQR KEERQAYLTN EYMQRMASVR AVPNQRAPPS GYFMTAVPQT QNHAA YPPS QIARLRPSR WTAQGARPHP FQNKPSAIRP GAPRVPFSTM RPASSQVPRV MSTQRVANTS TQTVGPRPAA AAAAAATPAV RTVPRYKYAA GVRNPQQHRN AQPQVTMQQL AVHVQGGQETL TASRLASAPP QKQKQMLGER LFPLIQAMHP TLAGKITGML LEIDNSEL LY MLESPELRS KVDEAVAVLQ AHQAKEATQK AVNSATGVPT V |
| Source | E.coli |
| Target Names | PABPC3 |
| Protein Names | Recommended name: Polyadenylate-binding protein 3 Short name= PABP-3 Short name= Poly(A)-binding protein 3 Alternative name(s): Testis-specific poly(A)-binding protein |
| Expression Region | 1-631 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final |



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.