



Recombinant Human WD repeat-containing protein C2orf44 (C2orf44)

Product Code	CSB-MP875684HU
Abbreviation	C2orf44
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H6R7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MELGKGKLLR TGLNALHQAV HPIHGLAWTD GNQVVLTDLR LHSGEVKFGD SKVIGQFECV CGLSWAPPVA DDTPVLLAVQ HEKHVTWVQL CPSPMESSKW LTSQTCEIRG SLPILPQGCV WHPKCAITV LTAQDVSIFP NVHSDDSQVK ADINTQGRIH CACWTQDGLR LVVAVGSSLH SYIWDSAQKT LHRCSSCLVF DVDSHVCSIT ATVDSQVAIA TELPLDKICG LNASETFNIP PNSKDMTPYA LPVIGEVRSMDKEATDSETN SEVSVSSSYL EPLDLTHIHF NQHKSEGNSL ICLRKKDYLT GTGQDSSHLV LVTFKKAVTM TRKVTIPGIL VPDLIAFNLK AHVVAVASNT CNILIYSVI PSSVPNIQQI RLEENTERPKG ICFLTDQLLL ILVGKQKLTDTTFLPSSKSD QYAISLIVRE IMLEEPSIT SGESQTTYST FSAPLNKANR KKLIESLSPD FCHQNKGLLL TVNTSSQNGR PGRTLKIEIQ SPLSSICDGS IALDAEPVTQ PASLPRHSST PDHTSTLEPP RLPQRKNLQS EKETYQLSKE VEILSRNLVE MQRCLSELTN RLHNGKKSSS VYPLSQDLPY VHIIYQKPYL LGPVVEKRAV LLCDGKLRLLS TVQQTFGLSL IEMLHDSHWI LLSADSEGFI PLTFTATQEI IIRDGSLRS DVFRDSFSHS PGAVSSLKVF TGLAAPSLDT TGCCNHVDGM A
Source	Mammalian cell
Target Names	WDCP
Protein Names	Recommended name: WD repeat-containing protein C2orf44
Expression Region	1-721
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.